

Important tubes at a glance																					Maximum performance of the centrifuge
Volume in ml	0.2-2	8x0.2	0.5	15	50	100	250	1.6-7	4-10	15	50	250/500	230/290	650/750	1,000	2,000	MTP	500	-	-	
H x W x D in mm																180x150 x100	86x128 x15				
Ø x L in mm	6-11x38	-	11x36	17x100	34x100	40x115	65x115	13x75/100	15x75/100	17x120	29x115/115	60x192/96x147	61x122/62x137	97x139/152	96x178			quadruple	-	-	

CENTRIFUGE MODELS	Maximum number of tubes per centrifuge										Maximum number of tubes per centrifuge										max. RCF	max. RPM min ⁻¹
HANDZENTRIFUGE	-	-	-	4	-	-	-	4	4	4	-	-	-	-	-	-	-	-	-	-	1,298	3,000
EBA 200/EBA 200 S	-	-	8	8	-	-	-	8	- / 8	4	-	-	-	-	-	-	-	-	-	-	3,461/6,158	8,000/8,000
EBA 270	-	-	6	6	-	-	-	6	6	-	-	-	-	-	-	-	-	-	-	-	2,254	4,000
EBA 280/EBA 280 S	-	-	12	12	6	-	-	12	8 / 12	6	3 / -	-	-	-	-	-	-	-	-	-	4,148/5,071	6,000
MIKRO 185	24	-	12	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	17,008	13,300
MIKRO 200/200 R	30	4x8	15	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	21,382	15,000
MIKRO 220/220 R	48 (60)	6x8	15	12	6	-	-	12	8* / 12	6	3 / -	-	-	-	-	-	-	-	-	-	31,514	18,000
ROTOFIX 32 A	(36)	-	24	32	6	4 (44x100 mm)	-	40/32	28/32	32	8	-	-	-	-	-	-	6	12	4,226	6,000	
UNIVERSAL 320/320 R	30	24x8	15	32	6	4 (44x100 mm)	-	40/32	28/32	32	8	-	4 x 200 ml	-	-	-	10	-	6	12	21,382	15,000
ROTINA 380/380 R	30 (60)	24x8	15	52	12	4 (44x100 mm)	4	64/76	52	36	16	-	4	-	-	-	12	-	8	-	24,400	15,000
ROTINA 420/420 R	30 (160)	48x8	15	72	16	12	4	104	72	52	20	-	4	4x800 ml	-	-	16	-	8	-	24,400	15,000
ROTANTA 460/460 R ROTANTA 460 RC/460 RF	30 (240)	48x8	15	148	16	12	6	196	148	96	40/28	4	8	4	-	-	24	4x450 ml	18	-	24,400	15,000
ROTOFIX 46	-	-	-	48	8	4	4	-	-	28	8	-	4	-	-	-	-	-	-	-	3,095	4,000
ROTIXA 500 RS	(336)	48x8	-	168	24	16	6	200	168	112	40	4	6	4	4	-	24	4	12	-	18,098	11,500
ROTO SILENTA 630 RS	-	-	-	180	36	24	-	180	180	138	48/36	6	6	6	6	6	-	12	-	-	6,520	6,000
MIKRO 220 Robotic	4 / - 24 / 12	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	18,516	13,000
ROTINA 380/380 R/380 RC Robotic	on request							48	48	on request	16	on request					6	on request			4,696	5,100
ROTANTA 460 Robotic	on request							80	80	on request	24	on request					12	on request			6,446	6,200
ROTOLAVIT II	-	-	-	24x5 ml	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	1,438	3,500